

Anti-Acy1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q03154

Item Number: CM0014263



Introduction

Rabbit polyclonal antibody targeting human aminoacylase-1 (ACY1). Validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to mouse and rat. Core reagent for N-acyl-amino acid hydrolysis and renal metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ACY1 Rabbit pAb
Remarks/Alias	ACY-1; ACY1D; HEL-S-5; ACY1
Species	Human
Gene ID (Human)	95
Gene ID	95
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-408 of human ACY1 (NP_001185824.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03154
Protein Weight	46kDa
Shipping	Ice bag