


Anti-Adam10 Polyclonal Antibody For Wb, If/Icc, Fc, Elisa - O14672

Item Number: CM0029841



Introduction

ADAM10 rabbit polyclonal antibody for Western blot, immunofluorescence/ICC, flow cytometry, and ELISA. Recognizes human and rat ADAM10, a key sheddase involved in neurobiology, inflammation, and cancer.

[Learn More](#)

Parameter	Value
Product Name	ADAM10 Rabbit pAb
Remarks/Alias	RAK; kuz; AD10; AD18; MADM; CD156c; CDw156; HsT18717; ADAM10
Species	Human
GeneID (Human)	102
GeneID	102
Immunogen	Recombinant protein corresponding to a sequence within amino acids 214-672 of human ADAM10 (NP_001101.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human, Rat
SWISS	O14672
Protein Weight	84kDa
Shipping	Ice bag