

Anti-Adam11 Polyclonal Antibody For Wb, If/Icc, Elisa - 075078

Item Number: CM0015418

Introduction

Rabbit polyclonal antibody targeting human ADAM11 (disintegrin and metalloproteinase domain-containing protein 11), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. For neuroscience studies of synaptic function and brain development.

[Learn More](#)

Parameter	Value
Product Name	ADAM11 Rabbit pAb
Remarks/Alias	MDC; ADAM11
Species	Human
GeneID (Human)	4185
GeneID	4185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-430 of human ADAM11 (NP_002381.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075078
Protein Weight	83kDa
Shipping	Ice bag