

Anti-Adam12 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43184

Item Number: CM0029174

Introduction

Rabbit polyclonal antibody against ADAM12 (Meltrin-alpha), validated for WB, IHC-P, IF/ICC, ELISA, reacts with human, mouse, rat. Involved in skeletal muscle regeneration and cell fusion. Ideal for muscle biology research.

[Learn More](#)

Parameter	Value
Product Name	ADAM12 Rabbit pAb
Remarks/Alias	MCMP, MLTN, CAR10, MLTNA, MCMPMitna, ADAM12-OT1, ADAM12
Species	Human
Gene ID	8038 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 589-738 of human ADAM12 (NP_067673.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43184
Protein Weight	100kDa
Shipping	Ice bag