

Anti-Adam15 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q13444

Item Number: CM0008672

Introduction

High-quality ADAM15 rabbit monoclonal antibody for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studies on cell adhesion, wound healing, and cancer metastasis.

[Learn More](#)

Parameter	Value
Product Name	ADAM15 Rabbit mAb
Remarks/Alias	MDC15; ADAM15
Species	Human
Gene ID (Human)	8751
Immunogen	Synthetic peptide corresponding to amino acids 200-300 of human ADAM15 (Q13444)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13444
Protein Weight	55kDa/88kDa/89kDa/90kDa/93kDa
Shipping	Ice bag