

Anti-Adam15 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13444

Item Number: CM0029727

Introduction

Rabbit polyclonal ADAM15 antibody for WB, IF/ICC, and ELISA applications, reactive with human, mouse, and rat samples. Suitable for wound healing, cell adhesion, and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	ADAM15 Rabbit pAb
Remarks/Alias	MDC15; ADAM15
Species	Human
GeneID (Human)	8751
GeneID	8751
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 540-700 of human ADAM15 (NP_003806.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13444
Protein Weight	93kDa
Shipping	Ice bag