

Anti-Adam23 Polyclonal Antibody For Wb, Elisa - 075077

Item Number: CM0014826

Introduction

Rabbit polyclonal ADAM23 antibody validated for WB and ELISA; reacts with human, mouse, and rat. Targets the non-catalytic metalloprotease-like protein involved in cell adhesion, highly expressed in brain.

[Learn More](#)

Parameter	Value
Product Name	ADAM23 Rabbit pAb
Remarks/Alias	MDC3; MDC-3; ADAM23
Species	Human
Gene ID (Human)	8745
Gene ID	8745
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-500 of human ADAM23 (NP_003803.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075077
Protein Weight	92kDa
Shipping	Ice bag