

## Anti-Adam28 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ukq2

Item Number: CM0014538

### Introduction

Rabbit polyclonal ADAM28 antibody validated for Western blot and ELISA. Reacts with human and mouse ADAM28, a disintegrin metalloproteinase involved in lymphocyte migration and ectodomain shedding.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | ADAM28 Rabbit pAb                                                                                                   |
| Remarks/Alias    | MDCL; MDC-L; eMDCII; ADAM 28; eMDC II; ADAM28                                                                       |
| Species          | Human                                                                                                               |
| GeneID (human)   | 10863                                                                                                               |
| GeneID           | 10863                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 199-440 of human ADAM28 (NP_055080.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse                                                                                                        |
| SWISS            | Q9UKQ2                                                                                                              |
| Protein Weight   | 87kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |