


Anti-Adam33 Polyclonal Antibody For Wb, Elisa - Q9Bz11

Item Number: CM0026267



Introduction

Rabbit polyclonal antibody against human ADAM33, validated for WB and ELISA, with cross-reactivity to mouse and rat. Targets the metalloprotease involved in angiogenesis, cell proliferation, and cytokine shedding.

[Learn More](#)

Parameter	Value
Product Name	ADAM33 Rabbit pAb
Remarks/Alias	C20orf153; DJ964F7.1; ADAM33
Species	Human
GeneID (Human)	80332
GeneID	80332
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human ADAM33 (NP_001269376.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZ11
Protein Weight	88kDa
Shipping	Ice bag