

Anti-Adam5 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Nvv9

Item Number: CM0025792

Introduction

Rabbit polyclonal antibody against ADAM5, a putative non-catalytic metalloprotease highly expressed in testis, validated for Western blot and ELISA in human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	ADAM5 Rabbit pAb
Remarks/Alias	ADAM5P; TMDCI; ADAM5
Species	Human
Gene ID (Human)	255926
Gene ID	255926
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 250-350 of human ADAM5 (Q6NVV9).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6NVV9
Protein Weight	47kDa
Shipping	Ice bag