

Anti-Adam8 Rabbit Polyclonal Antibody For Wb And Elisa - P78325

Item Number: CM0016886



Introduction

Rabbit polyclonal antibody targeting human ADAM8, suitable for Western blot and ELISA. Cross-reacts with mouse and rat ADAM8. Ideal for investigating leukocyte extravasation, metalloprotease activity, and inflammatory processes.

[Learn More](#)

Parameter	Value
Product Name	ADAM8 Rabbit pAb
Remarks/Alias	MS2; CD156; CD156a; ADAM8
Species	Human
Gene ID (Human)	101
Gene ID	101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-300 of human ADAM8 (NP_001100.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78325
Protein Weight	89kDa
Shipping	Ice bag