

Anti-Adam9 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q13443

Item Number: CM0007495

Introduction

Rabbit monoclonal ADAM9 antibody validated for WB, IF-P, IHC-P, ELISA; reacts with human, mouse, monkey. ADAM9 is a metalloprotease involved in cell adhesion, migration, and tumor angiogenesis. Ideal for cancer and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	ADAM9 Rabbit mAb
Remarks/Alias	MCMP; MDC9; CORD9; Mltng; ADAM9
Species	Human
GeneID (human)	8754
GeneID	8754
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-210 of human ADAM9 (NP_003807.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Monkey
SWISS	Q13443
Protein Weight	91kDa
Shipping	Ice bag