

Anti-Adamts1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uhi8

Item Number: CM0016628



Introduction

Rabbit polyclonal antibody against human ADAMTS1 (METH1). Validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Molecular weight: 105kDa. Suitable for extracellular matrix and angiogenesis research.

[Learn More](#)

Parameter	Details
Product Name	ADAMTS1 Rabbit pAb
Remarks / Alias	C3-C5; METH1; ADAMTS1
Species	Human
GeneID (Human)	9510
GeneID	9510
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 718-967 of human ADAMTS1 (NP_008919.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHI8
Protein Weight	105kDa
Shipping	Ice bag