

Anti-Adamts13 Polyclonal Antibody For Wb, If/Icc, Elisa - Q76Lx8

Item Number: CM0026274



Introduction

Rabbit polyclonal antibody against human ADAMTS13 (vWF-cleaving protease). Suitable for WB, IF/ICC, ELISA. Detects human, mouse, and rat ADAMTS13. Research reagent for hemostasis and thrombotic disease studies.

[Learn More](#)

Parameter	Value
Product Name	ADAMTS13 Rabbit pAb
Remarks/Alias	VWFCP; C9orf8; vWF-CP; ADAM-TS13; ADAMTS-13; ADAMTS13
Species	Human
Gene ID	11093
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-300 of human ADAMTS13 (NP_620594.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q76LX8
Protein Weight	154kDa
Shipping	Ice bag