

Anti-Adamts2 Rabbit Polyclonal Antibody For Wb, Elisa - 095450

Item Number: CM0017949

Introduction

High-quality anti-ADAMTS2 rabbit polyclonal antibody for Western blot and ELISA. Reacts with human and rat ADAMTS2, a procollagen N-proteinase essential for collagen fibril assembly. Suitable for extracellular matrix and connective tissue research.

[Learn More](#)

Parameter	Value
Product Name	ADAMTS2 Rabbit pAb
Remarks/Alias	NPI; PNPI; PCINP; PCPNI; PCI-NP; PC I-NP; ADAM-TS2; ADAMTS-2; ADAMTS-3; EDSDERMS; ADAMTS2
Species	Human
GeneID (Human)	9509
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1029-1213 of human ADAMTS2 (095450)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	095450
Protein Weight	135kDa
Shipping	Ice bag