


Anti-Adamts4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 075173

Item Number: CM0024311



Introduction

High-quality rabbit polyclonal antibody targeting human ADAMTS4 (aggrecanase-1). Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for research on arthritis, cartilage degradation, and Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	ADAMTS4 Rabbit pAb
Remarks/Alias	ADMP-1; ADAMTS-2; ADAMTS-4; ADAMTS4
Species	Human
Gene ID (Human)	9507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 74-318 of human ADAMTS4 (NP_005090.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075173
Protein Weight	90kDa
Shipping	Ice bag