


Anti-Adamts5 Polyclonal Antibody For Wb, Elisa - Q9Una0

Item Number: CM0024878



Introduction

Rabbit polyclonal antibody targeting human ADAMTS5 (aggrecanase-2) for Western blot (WB) and ELISA applications. Recognizes human, mouse, and rat ADAMTS5, with a molecular weight of approximately 101 kDa.

[Learn More](#)

Parameter	Value
Product Name	ADAMTS5 Rabbit pAb
Remarks/Alias	ADAMTS5; ADAM-TS 11; ADAM-TS 5; ADAM-TS5; ADAMTS-11; ADAMTS-5; ADAMTS11; ADMP-2; ADAM-TS11
Species	Human
GeneID (Human)	11096
GeneID	11096
Immunogen	Synthetic peptide corresponding to amino acids 601-700 of human ADAMTS5 (NP_008969.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNA0
Protein Weight	101kDa
Shipping	Ice bag