


Anti-Adamts5 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Una0

Item Number: CM0005631



Introduction

This recombinant rabbit monoclonal antibody targets human ADAMTS5 (aggrecanase-2) with cross-reactivity to mouse, validated for WB and ELISA. Ideal for studies on extracellular matrix degradation, tissue remodeling, and protease function.

[Learn More](#)

Parameter	Value
Product Name	ADAMTS5 Rabbit mAb
Remarks/Alias	ADMP-2; ADAM-TS5; ADAMTS-5; ADAMTS11; ADAM-TS 5; ADAMTS-11; ADAM-TS 11; ADAMTS5
Species	Human
GeneID (Human)	11096
GeneID	11096
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 601-700 of human ADAMTS5 (NP_008969.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UNA0
Protein Weight	102 kDa
Shipping	Ice bag