

Anti-Adamts7 Rabbit Polyclonal Antibody For Wb & Elisa - Q9Ukp4

Item Number: CM0019040

Introduction

High-quality rabbit polyclonal antibody against human ADAMTS7, validated for WB and ELISA. Suitable for studying metalloprotease function in extracellular matrix remodeling and related diseases.

[Learn More](#)

Parameter	Value
Product Name	ADAMTS7 Rabbit pAb
Remarks/Alias	ADAM-TS7; ADAMTS-7; ADAM-TS 7; ADAMTS7
Species	Human
GeneID (Human)	11173
GeneID	11173
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1600-1686 of human ADAMTS7 (NP_055087.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UKP4
Protein Weight	184kDa
Shipping	Ice bag