

Anti-Adamtsl1 Rabbit Polyclonal Antibody For Wb And Elisa - Q8N6G6

Item Number: CM0028517

Introduction

High-quality rabbit polyclonal antibody targeting human ADAMTSL1, suitable for WB and ELISA. Cross-reacts with mouse and rat. Recombinant protein immunogen. For skeletal muscle and extracellular matrix studies.

[Learn More](#)

Parameter	Value
Product Name	ADAMTSL1 Rabbit pAb
Remarks/Alias	C9orf94; PUNCTIN; ADAMTSR1; ADAMTSL-1; ADAMTSL1
Species	Human
GeneID (Human)	92949
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1483-1762 of human ADAMTSL1 (NP_001035362.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N6G6
Protein Weight	62kDa
Shipping	Ice bag