


Anti-Adar1 Rabbit Monoclonal Antibody For Wb, Elisa - P55265

Item Number: CM0002402



Introduction

High-quality anti-ADAR1 rabbit monoclonal antibody validated for Western blot and ELISA applications. Specific to human ADAR1 (P55265), a 136kDa double-stranded RNA-specific adenosine deaminase involved in A-to-I RNA editing.

[Learn More](#)

Parameter	Value
Product Name	ADAR1 Rabbit mAb
Remarks/Alias	DSH; AGS6; G1P1; IFI4; P136; ADAR1; DRADA; DSRAD; IFI-4; K88DSRBP
Species	Human
GeneID (Human)	103
GeneID	103
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human ADAR1 (P55265).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P55265
Protein Weight	136kDa
Shipping	Ice bag