


Anti-Adar1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P55265

Item Number: CM0027006



Introduction

High-quality rabbit polyclonal antibody targeting human ADAR1, validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for investigating A-to-I RNA editing, antiviral immunity, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	ADAR1 Rabbit pAb
Remarks/Alias	DSH; AGS6; G1P1; IFI4; P136; ADAR1; DRADA; DSRAD; IFI-4; K88DSRBP
Species	Human
GeneID (human)	103
GeneID	103
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-299 of human ADAR1 (NP_001102.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55265
Protein Weight	136kDa
Shipping	Ice bag