

Anti-Adarb1 Rabbit Polyclonal Antibody For Wb, Elisa - P78563

Item Number: CM0020514



Introduction

High-quality anti-ADARB1 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat ADARB1 (81 kDa). Ideal for A-to-I RNA editing and neurological disease research.

[Learn More](#)

Parameter	Value
Product Name	ADARB1 Rabbit pAb
Remarks/Alias	RED1, ADAR2, DRABA2, DRADA2, NEDHYMS, ADARB1
Species	Human
GeneID (Human)	104
GeneID	104
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-380 of human ADARB1 (NP_001103.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78563
Protein Weight	81kDa
Shipping	Ice bag