


Anti-Adcy10/Sac Rabbit Polyclonal Antibody For Wb And Elisa - Q96Pn6

Item Number: CM0019090



Introduction

Rabbit polyclonal antibody targeting human ADCY10 (soluble adenylyl cyclase), validated for WB and ELISA applications. Reacts with human and mouse samples; ideal for cAMP signaling and fertility studies.

[Learn More](#)

Parameter	Value
Product Name	ADCY10 Rabbit pAb
Remarks/Alias	SAC; HCA2; SACI; Sacy; hsAC; HEL-S-7a; ADCY10
Species	Human
GeneID (Human)	55811
GeneID	55811
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human ADCY10 (NP_060887.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96PN6
Protein Weight	187kDa
Shipping	Ice bag