

Anti-Adcy2 Rabbit Polyclonal Antibody For Wb, Elisa - Q08462

Item Number: CM0029924

Introduction

Rabbit polyclonal antibody against human ADCY2, validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for studying cAMP signaling, GPCR pathways, and adrenal/brain tissue expression.

[Learn More](#)

Parameter	Value
Product Name	ADCY2 Rabbit pAb
Remarks/Alias	AC2; HBAC2; ADCY2
Species	Human
Gene ID (Human)	108
Gene ID	108
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 430-600 of human ADCY2 (NP_065433.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08462
Protein Weight	124kDa
Shipping	Ice bag