

Anti-Adcyap1R1 (Pac1) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P41586

Item Number: CM0025214



Introduction

Rabbit polyclonal antibody against human ADCYAP1R1 (PAC1/PACAP receptor), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse, rat. Ideal for GPCR signaling and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	ADCYAP1R1 Rabbit pAb
Remarks/Alias	PAC1; PAC1R; PACAPR; PACAPRI; ADCYAP1R1
Species	Human
GeneID (Human)	117
GeneID	117
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-140 of human ADCYAP1R1 (NP_001109.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41586
Protein Weight	53kDa
Shipping	Ice bag