


Anti-Adenosine Deaminase (Ada) Rabbit Monoclonal Antibody For Wb, Elisa - P00813

Item Number: CM0008309



Introduction

Adenosine deaminase (ADA) rabbit monoclonal antibody validated for WB and ELISA, showing cross-reactivity with human, mouse, and rat ADA. Ideal for purine metabolism, immunology, and disease research.

[Learn More](#)

Parameter	Value
Product Name	Adenosine Deaminase (ADA) Rabbit mAb
Remarks/Alias	ADA1; Adenosine Deaminase (ADA)
Species	Human
GeneID (Human)	100
GeneID	100
Immunogen	Synthetic Peptide corresponding to amino acids 250-350 of human ADA (NP_000013.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00813
Protein Weight	41kDa
Shipping	Ice bag