


Anti-Adenosine Deaminase (Ada) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P00813

Item Number: CM0013439



Introduction

Rabbit polyclonal antibody targeting human Adenosine Deaminase (ADA). Validated for WB, IF/ICC, ELISA. Cross-reacts with human and rat. Suitable for purine metabolism, T-cell activation, and adenosine signaling research.

[Learn More](#)

Parameter	Value
Product Name	Adenosine Deaminase (ADA) Rabbit pAb
Remarks/Alias	ADA1; Adenosine Deaminase (ADA)
Species	Human
Gene ID (Human)	100
Gene ID	100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-363 of human Adenosine Deaminase (ADA) (NP_000013.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P00813
Protein Weight	41kDa
Shipping	Ice bag