


Anti-Adenosylhomocysteinase (Ahcy/Sahh) Rabbit Monoclonal Antibody For Wb And Elisa - P23526

Item Number: CM0011887



Introduction

High-quality rabbit monoclonal antibody against human AHCY (SAHH), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying one-carbon metabolism and homocysteine pathways.

[Learn More](#)

Product Name	AHCY Rabbit mAb
Remarks/Alias	SAHH; adoHcyase; AHCY
Species	Human
GeneID (Human)	191
GeneID	191
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 300-400 of human AHCY (P23526).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23526
Protein Weight	48kDa
Shipping	Ice bag