

Anti-Adenylate Kinase 4 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P27144

Item Number: CM0004223



Introduction

Rabbit monoclonal antibody targeting human Adenylate kinase 4 (AK4), suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with rat. Ideal for studies on mitochondrial nucleotide homeostasis, AMPK signaling, and oxidative stress.

[Learn More](#)

Parameter	Value
Product Name	Adenylate kinase 4 Rabbit mAb
Remarks/Alias	AK3; AK 4; AK3L1; AK3L2; Adenylate kinase 4
Species	Human
Gene ID (Human)	205
Gene ID	205
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human Adenylate kinase 4 (P27144)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P27144
Protein Weight	25kDa
Shipping	Ice bag