

Anti-Adgra3 Rabbit Polyclonal Antibody For Wb, Elisa - Q8lwk6

Item Number: CM0017170

Introduction

Rabbit polyclonal antibody against human ADGRA3 (GPR125) for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for studying adhesion GPCRs and planar cell polarity pathways.

[Learn More](#)

Parameter	Value
Product Name	ADGRA3 Rabbit pAb
Remarks/Alias	PGR21; TEM5L; GPR125; ADGRA3
Species	Human
Gene ID (Human)	166647
Gene ID	166647
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 370-600 of human ADGRA3 (NP_660333.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8lwk6
Protein Weight	146kDa
Shipping	Ice bag