


Anti-Adgrl4 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hbw9

Item Number: CM0014869



Introduction

Rabbit polyclonal antibody targeting human ADGRL4 (ELTD1/ETL) for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Validated for angiogenesis research, 78 kDa, Swiss-Prot Q9HBW9.

[Learn More](#)

Parameter	Value
Product Name	ADGRL4 Rabbit pAb
Remarks/Alias	ETL; ELTD1; KPG_003; ADGRL4
Species	Human
GeneID (Human)	64123
GeneID	64123
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 110-406 of human ADGRL4 (NP_071442.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HBW9
Protein Weight	78kDa
Shipping	Ice bag