

Anti-Adh1A/Adh1B/Adh1C Rabbit Polyclonal Antibody For Wb - P07327 / P00325 / P00326

Item Number: CM0021632



Introduction

CamelBio rabbit polyclonal antibody targeting ADH1A/ADH1B/ADH1C alcohol dehydrogenases. Validated for WB, IF-P, ELISA in human, mouse, rat. Ideal for alcohol metabolism and liver function studies.

[Learn More](#)

Parameter	Value
Product Name	ADH1A/ADH1B/ADH1C Rabbit pAb
Remarks/Alias	ADH1A/ADH1B/ADH1C
Species	Human
Gene ID (Human)	124, 125, 126
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07327, P00325, P00326
Shipping	Ice bag