

Anti-Adh4 Rabbit Polyclonal Antibody For Wb, Elisa - P08319

Item Number: CM0028576



Introduction

CamelBio offers Anti-ADH4 Rabbit Polyclonal Antibody validated for WB and ELISA in human samples, with mouse cross-reactivity. Ideal for retinol metabolism and alcohol dehydrogenase research. Explore reliable IVD raw materials.

[Learn More](#)

Parameter	Value
Product Name	ADH4 Rabbit pAb
Remarks/Alias	ADH-2; HEL-S-4; ADH4
Species	Human
GeneID (Human)	127
GeneID	127
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-380 of human ADH4 (NP_000661.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P08319
Protein Weight	40kDa
Shipping	Ice bag