

## Anti-Adh5/Gsnor Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P11766

Item Number: CM0007255



### Introduction

Rabbit monoclonal antibody targeting ADH5/GSNOR (SWISS P11766), validated for WB, IF-P, ELISA. Detects human, mouse, rat ADH5/GSNOR. For alcohol metabolism, formaldehyde detoxification, and S-nitrosylation studies.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ADH5/GSNOR Rabbit mAb                                                                                                  |
| Remarks/Alias    | FDH; ADHX; ADH-3; AMEDS; BMFS7; FALDH; GSNOR; GSH-FDH; HEL-S-60p; ADH5/GSNOR                                           |
| Species          | Human                                                                                                                  |
| Gene ID (Human)  | 128                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-374 of human ADH5/GSNOR (NP_000662.3). |
| Source           | Rabbit                                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                                  |
| Application      | WB, IF-P, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                      |
| SWISS            | P11766                                                                                                                 |
| Protein Weight   | 40kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                |