

Anti-Adhfe1 Polyclonal Antibody For Wb, Elisa - Q8lww8

Item Number: CM0022711

Introduction

Rabbit polyclonal antibody targeting human ADHFE1 (HOT/ADH8). Applications: WB, ELISA; cross-reacts with mouse. Recombinant immunogen (aa 41-255). Ideal for studying mitochondrial GHB metabolism. Molecular weight: 50 kDa.

[Learn More](#)

Parameter	Specification
Product Name	ADHFE1 Rabbit pAb
Remarks/Alias	HOT; ADH8; HMFT2263
Species	Human
Gene ID (Human)	137872
Gene ID	137872
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 41-255 of human ADHFE1 (NP_653251.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8IWW8
Protein Weight	50kDa
Shipping	Ice bag