


Anti-Adiponectin Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15848

Item Number: CM0027876



Introduction

Validated rabbit polyclonal antibody against human adiponectin (SWISS Q15848) for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for metabolic disease research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Adiponectin Rabbit pAb
Remarks/Alias	ACDC; ADPN; APM1; APM-1; GBP28; ACRP30; ADIPQTL1; Adiponectin
Species	Human
Gene ID	9370
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-244 of human Adiponectin (NP_004788.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15848
Protein Weight	26 kDa
Shipping	Ice bag