

Anti-Adiponectin Rabbit Monoclonal Antibody For Wb And Elisa - Q15848

Item Number: CM0007337



Introduction

Adiponectin Rabbit mAb (Q15848) validated for Western blotting and ELISA, with cross-reactivity to mouse. A key adipokine in metabolism and insulin sensitivity, ideal for diabetes and obesity research.

[Learn More](#)

Parameter	Value
Product Name	Adiponectin Rabbit mAb
Remarks/Alias	ACDC; ADPN; APM1; APM-1; GBP28; ACRP30; ADIPQTL1; Adiponectin
Species	Human
Gene ID (Human)	9370
Gene ID	9370
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-244 of human Adiponectin. (NP_004788.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q15848
Protein Weight	26kDa
Shipping	Ice bag