


Anti-Adipor1 Polyclonal Antibody For Wb, Elisa - Q96A54

Item Number: CM0016120



Introduction

Rabbit polyclonal antibody targeting human ADIPOR1 (adiponectin receptor 1, 43 kDa). Validated for WB and ELISA with cross-reactivity in mouse and rat. Ideal for metabolic disease research, diabetes, and obesity studies.

[Learn More](#)

Parameter	Value
Product Name	ADIPOR1 Rabbit pAb
Remarks/Alias	CGI45; PAQR1; ACDCR1; CGI-45; TESBP1A; ADIPOR1
Species	Human
Gene ID (Human)	51094
Gene ID	51094
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-59 of human Adiponectin Receptor 1 (NP_057083.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96A54
Protein Weight	43kDa
Shipping	Ice bag