


Anti-Adipor1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96A54

Item Number: CM0024169



Introduction

Rabbit polyclonal antibody targeting ADIPOR1, validated for WB and ELISA with cross-reactivity to Mouse and Rat. Ideal for studying adiponectin signaling and metabolic regulation in glucose and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ADIPOR1 Rabbit pAb
Remarks/Alias	CGI45; PAQR1; ACDCR1; CGI-45; TESBP1A; ADIPOR1
Species	Human
Gene ID (Human)	51094
Gene ID	51094
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96A54
Protein Weight	43kDa
Shipping	Ice bag