

Anti-Adipor2 Polyclonal Antibody For Wb, Elisa - Q86V24

Item Number: CM0017076



Introduction

ADIPOR2 Rabbit pAb detects human, mouse, rat adiponectin receptor 2 (44kDa) via WB and ELISA. Ideal for metabolic studies on glucose and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	ADIPOR2 Rabbit pAb
Remarks/Alias	PAQR2; ACDCR2; ADIPOR2
Species	Human
GeneID (Human)	79602
GeneID	79602
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human ADIPOR2 (NP_078827.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86V24
Protein Weight	44kDa
Shipping	Ice bag