

Anti-Adk Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P55263

Item Number: CM0013536



Introduction

High-quality rabbit polyclonal antibody against human adenosine kinase (ADK). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Ideal for adenosine metabolism and RNA modification research.

[Learn More](#)

Parameter	Value
Product Name	ADK Rabbit pAb
Remarks/Alias	AK; ADK
Species	Human
GeneID (Human)	132
GeneID	132
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 200-345 of human ADK (NP_001114.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55263
Protein Weight	41kDa
Shipping	Ice bag