

Anti-Adm Rabbit Polyclonal Antibody For Wb And Elisa - P35318

Item Number: CM0026462



Introduction

ADM Rabbit pAb is a rabbit polyclonal antibody raised against human ADM peptide. Validated for Western blot and ELISA, it detects human ADM and cross-reacts with mouse ADM. Suitable for vasodilation, blood pressure, and fluid balance research.

[Learn More](#)

Parameter	Specification
Product Name	ADM Rabbit pAb
Remarks/Alias	AM; PAMP; ADM
Species	Human
Gene ID (Human)	133
Gene ID	133
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human ADM (NP_001115.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P35318
Protein Weight	20kDa
Shipping	Ice bag