

Anti-Adnp Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H2P0

Item Number: CM0027333

Introduction

Rabbit polyclonal antibody against human ADNP (Q9H2P0) for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse. Targets activity-dependent neuroprotector homeobox protein involved in transcriptional regulation and WNT signaling.

[Learn More](#)

Parameter	Value
Product Name	ADNP Rabbit pAb
Remarks/Alias	ADNP1; HVDAS; MRD28; ADNP
Species	Human
GeneID (Human)	23394
GeneID	23394
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 843-1102 of human ADNP (NP_056154.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H2P0
Protein Weight	124kDa
Shipping	Ice bag