


Anti-Adora1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P30542

Item Number: CM0004481



Introduction

Anti-ADORA1 rabbit mAb for WB, IHC-P, ELISA in human, mouse, rat. Targets adenosine receptor A1, a GPCR mediating inhibitory signaling. Research use only.

[Learn More](#)

Parameter	Value
Product Name	ADORA1 Rabbit mAb
Remarks/Alias	RDC7; ADORA1
Species	Human
Gene ID (Human)	134
Gene ID	134
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human ADORA1 (P30542)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30542
Protein Weight	37kDa
Shipping	Ice bag