


## Anti-Adora1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P30542

Item Number: CM0020142



### Introduction

Rabbit polyclonal antibody targeting human ADORA1 (adenosine receptor A1), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Detects a band of 37 kDa. Ideal for GPCR signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                                       |
|------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ADORA1 Rabbit pAb                                                                                                           |
| Remarks/Alias    | RDC7; ADORA1                                                                                                                |
| Species          | Human                                                                                                                       |
| Gene ID (Human)  | 134                                                                                                                         |
| Gene ID          | 134                                                                                                                         |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 250-326 of human ADORA1 (NP_000665.1). |
| Source           | Rabbit                                                                                                                      |
| Category         | Polyclonal Antibodies                                                                                                       |
| Application      | WB, IHC-P, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                           |
| SWISS            | P30542                                                                                                                      |
| Protein Weight   | 37kDa                                                                                                                       |
| Shipping         | Ice bag                                                                                                                     |