


Anti-Adora2B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P29275

Item Number: CM0024383



Introduction

Rabbit polyclonal antibody targeting human adenosine receptor A2b (ADORA2B), validated for western blot, immunofluorescence/ICC, and ELISA, with cross-reactivity to mouse and rat. Ideal for adenosine signaling and GPCR research.

[Learn More](#)

Parameter	Value
Product Name	ADORA2B Rabbit pAb
Remarks/Alias	ADORA2; ADORA2B
Species	Human
GeneID (Human)	136
GeneID	136
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human ADORA2B (NP_000667.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P29275
Protein Weight	36kDa
Shipping	Ice bag