

Anti-Adora3 Polyclonal Antibody For Wb, Elisa - P0Dms8

Item Number: CM0022160



Introduction

Rabbit polyclonal antibody targeting human ADORA3 (A3 adenosine receptor). Validated for Western blot (WB) and ELISA applications. Cross-reacts with mouse and rat. Suitable for GPCR signaling and immunology research.

[Learn More](#)

Parameter	Value
Product Name	ADORA3 Rabbit pAb
Remarks/Alias	A3AR; ADORA3
Species	Human
GeneID (Human)	140
GeneID	140
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human ADORA3 (NP_000668.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P0DMS8
Protein Weight	36kDa
Shipping	Ice bag