

Anti-Adra1A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P35348

Item Number: CM0014233



Introduction

Rabbit polyclonal antibody targeting human ADRA1A (Alpha-1A adrenergic receptor), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Useful for studying catecholamine signaling and cardiac function.

[Learn More](#)

Parameter	Value
Product Name	ADRA1A Rabbit pAb
Remarks/Alias	ADRA1C; ADRA1L1; ALPHA1AAR; ADRA1A
Species	Human
GeneID (human)	148
GeneID	148
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 346-475 of human ADRA1A (NP_150646.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35348
Protein Weight	51kDa
Shipping	Ice bag