

Anti-Adra1A Rabbit Monoclonal Antibody For Wb, Elisa - P35348

Item Number: CM0009539



Introduction

Rabbit monoclonal antibody targeting human ADRA1A (Alpha-1A adrenergic receptor), suitable for Western blotting (WB) and ELISA. Detects ADRA1A at ~51 kDa in human and mouse samples. Ideal for GPCR signaling research.

[Learn More](#)

Parameter	Value
Product Name	ADRA1A Rabbit mAb
Remarks/Alias	ADRA1C; ADRA1L1; ALPHA1AAR; ADRA1A
Species	Human
GeneID	148
Immunogen	Synthetic peptide corresponding to amino acids 150-250 of human ADRA1A (P35348)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35348
Protein Weight	51kDa
Shipping	Ice bag